

Phospho-Src Family-Y416 Rabbit Pab For Wb And Elisa - P12931 / P06241 / P06239 / P07947 / P08631 / P07948

Item Number: CM0000909



Introduction

Rabbit polyclonal antibody targeting phosphorylated Tyr416 of SRC family kinases (SRC, FYN, YES1, LYN, HCK, LCK). Validated for WB and ELISA; cross-reacts with human, mouse, and rat. Ideal for studying cell signaling in cancer and immunity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-SRC Family-Y416 Rabbit pAb
Remarks/Alias	ASV; SRC1; THC6; c-SRC; p60-Src; SRC; Phospho-SRC Family-Y416
Species	Human
GeneID (Human)	6714/2534/3932/7525/3055/4067
GeneID	6714/2534/3932/7525/3055/4067
Immunogen	A synthetic phosphorylated peptide around Y416 of human SRC Family (NP_005408.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12931 / P06241 / P06239 / P07947 / P08631 / P07948
Protein Weight	60kDa
Shipping	Ice bag