


Phospho-Src-Y419 Rabbit Monoclonal Antibody For Wb And Elisa - P12931

Item Number: CM0000997



Introduction

Phospho-specific rabbit monoclonal antibody targeting Src Tyr419. Validated for Western blot and ELISA against human and mouse samples. Ideal for investigating Src kinase activation and signal transduction.

[Learn More](#)

Parameter	Specification
Product Name	Phospho-Src-Y419 Rabbit mAb
Remarks/Alias	ASV, SRC1, THC6, c-SRC, p60-Src, Phospho-Src-Y419
Species	Human
Gene ID (Human)	6714
Gene ID	6714
Immunogen	A synthetic phosphorylated peptide around Y419 of human SRC (P12931)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P12931
Protein Weight	60kDa
Shipping	Ice bag