


Phospho-Src-Y529 Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P12931

Item Number: CM0001685



Introduction

Phospho-specific rabbit polyclonal antibody targeting Src at Tyr529, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, rat. Ideal for studying Src kinase regulation and signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Src-Y529 Rabbit pAb
Remarks/Alias	ASV; SRC1; THC6; c-SRC; p60-Src; Phospho-Src-Y529
Species	Human
GeneID (Human)	6714
GeneID	6714
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12931
Protein Weight	60kDa
Shipping	Ice bag