

## Phospho-Stat3-S727 Rabbit Polyclonal Antibody For Wb, Ip, Ihc-P, Elisa - P40763

Item Number: CM0000534



### Introduction

Phospho-specific rabbit polyclonal antibody recognizing STAT3 phosphorylated at Ser727. Validated for Western blot, IP, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for studying STAT3 signaling and transcription regulation.

[Learn More](#)

| Parameter        | Details                                                                                                     |
|------------------|-------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-STAT3-S727 Rabbit pAb                                                                               |
| Remarks/Alias    | APRF; HIES; ADMIO; ADMIO1; Phospho-STAT3-S727                                                               |
| Species          | Human                                                                                                       |
| Gene ID (Human)  | 6774                                                                                                        |
| Gene ID          | 6774                                                                                                        |
| Immunogen        | Synthetic peptide: A synthetic phosphorylated peptide around S727 of human Phospho-STAT3-S727 (NP_644805.1) |
| Source           | Rabbit                                                                                                      |
| Category         | Phosphorylated Antibodies                                                                                   |
| Application      | WB, IP, IHC-P, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                           |
| SWISS            | P40763                                                                                                      |
| Protein Weight   | 88 kDa                                                                                                      |
| Shipping         | Ice bag                                                                                                     |