

Phospho-Stat3-Y705 Rabbit Mab For Wb, If-P, Ihc-P, Elisa - P40763

Item Number: CM0000686



Introduction

High-specificity Phospho-STAT3-Y705 Rabbit Monoclonal Antibody for WB, IF-P, IHC-P, and ELISA. Validated in human, mouse, and rat. Ideal for studying STAT3 signaling and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT3-Y705 Rabbit mAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1; Phospho-STAT3-Y705
Species	Human
GeneID (Human)	6774
GeneID	6774
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around Y705 of human STAT3 (NP_644805.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88kDa
Shipping	Ice bag