

Phospho-Stat3-Y705 Rabbit Mab - P40763

Item Number: CM0000979



Introduction

High-quality phospho-STAT3 (Tyr705) rabbit monoclonal antibody validated for WB, IP, ELISA. Specific for human STAT3, detects phosphorylation at Y705 activation site. Suitable for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT3-Y705 Rabbit mAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1; Phospho-STAT3-Y705
Species	Human
Gene ID	6774
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y705 of human STAT3 (NP_644805.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P40763
Protein Weight	88 kDa
Shipping	Ice bag