

Phospho-Stat3-Y705 Rabbit Pab For Wb, Ihc-P, Elisa - P40763

Item Number: CM0001549



Introduction

Phospho-specific rabbit polyclonal antibody targeting STAT3 phosphorylated at Y705, validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for studying JAK/STAT signaling and immune response pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT3-Y705 Rabbit pAb
Remarks/Alias	APRF; HIES; ADMIO; ADMIO1; Phospho-STAT3-Y705
Species	Human
GeneID (Human)	6774
GeneID	6774
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40763
Protein Weight	88kDa
Shipping	Ice bag