


Phospho-Stat5A-Y694 Rabbit Mab For Wb, Ihc-P, Elisa - P42229

Item Number: CM0001136



Introduction

Rabbit monoclonal antibody against phosphorylated STAT5A at Tyr694. Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for studying STAT5A signaling in cytokine responses and lactation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT5A-Y694 Rabbit mAb
Remarks/Alias	MGF; STAT5; Phospho-STAT5A-Y694
Species	Human
GeneID	6776
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around Y694 of human Stat5 (P42229).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42229
Protein Weight	91kDa
Shipping	Ice bag