

Phospho-Stat6-Y641 Rabbit Pab - P42226

Item Number: CM0001747



Introduction

Phospho-STAT6-Y641 Rabbit pAb specifically recognizes STAT6 phosphorylated at Y641. Suitable for Western blot, IHC-P, and ELISA in human, mouse, rat. An essential tool for investigating IL-4/STAT6 signaling pathways and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STAT6-Y641 Rabbit pAb
Remarks/Alias	STAT6B; STAT6C; D12S1644; IL-4-STAT; Phospho-STAT6-Y641
Species	Human
GeneID (human)	6778
GeneID	6778
Immunogen	A synthetic phosphorylated peptide around Y641 of human STAT6 (NP_003144.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42226
Protein Weight	94kDa
Shipping	Ice bag