

Phospho-Stathmin 1-S38 Rabbit Pab - P16949

Item Number: CM0001020

Introduction

Phospho-specific rabbit polyclonal antibody targeting Stathmin 1 pS38. Validated for WB and ELISA; reacts with human, mouse, and rat. Ideal for microtubule dynamics and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Stathmin 1-S38 Rabbit pAb
Remarks/Alias	Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215; Phospho-Stathmin 1-S38
Species	Human
Gene ID (Human)	3925
Gene ID	3925
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S38 of human Stathmin 1 (NP_001138926.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16949
Protein Weight	17kDa
Shipping	Ice bag