

Phospho-Sting/Tmem173-S365 Rabbit Mab For Wb, Elisa - Q3Tbt3

Item Number: CM0001556



Introduction

Phospho-STING/TMEM173-S365 Rabbit monoclonal antibody for Western blot and ELISA. Specifically detects mouse STING phosphorylated at serine 365, a critical activation marker in the cGAS-STING innate immune pathway.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STING/TMEM173-S365 Rabbit mAb
Remarks/Alias	ERIS; MPYS; Mita; STING; Tmem173; STING-beta; 2610307008Rik
Species	Mouse
GeneID (Human)	340061
GeneID	72512
Immunogen	Synthetic phosphorylated peptide around S365 of mouse STING (NP_082537.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q3TBT3
Protein Weight	43kDa
Shipping	Ice bag