

Phospho-Sting/Tmem173-S365 Rabbit Pab - Q3Tbt3

Item Number: CM0001312



Introduction

Phospho-STING/TMEM173-S365 Rabbit pAb detects phosphorylated mouse STING at S365 with cross-reactivity to human. Suitable for WB and ELISA, this antibody aids innate immunity and autophagy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STING/TMEM173-S365 Rabbit pAb
Remarks/Alias	ERIS; MPYS; Mita; STING; Tmem173; STING-beta; 2610307O08Rik; Phospho-STING-S365
Species	Mouse
GeneID (Human)	340061
GeneID	72512
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S365 of mouse STING (NP_082537.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q3TBT3
Protein Weight	43kDa
Shipping	Ice bag