

Phospho-Sting/Tmem173-S366 Rabbit Pab For Wb And Elisa - Q86Wv6

Item Number: CM0001798



Introduction

Phospho-STING/TMEM173-S366 rabbit polyclonal antibody for Western blotting and ELISA, targeting human STING/TMEM173 phosphorylated at Ser366, a pivotal adaptor in innate immune signaling and interferon induction.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STING/TMEM173-S366 Rabbit pAb
Remarks/Alias	ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; TMEM173; STING-beta; Phospho-STING/TMEM173-S366
Species	Human
Gene ID (Human)	340061
Gene ID	340061
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around S366 of human STING/TMEM173 (NP_938023.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q86WV6
Protein Weight	42kDa
Shipping	Ice bag