


Phospho-Stk11 (S428) Rabbit Polyclonal Antibody For Wb, Elisa - Q15831

Item Number: CM0001528



Introduction

Rabbit polyclonal antibody targeting phospho-STK11 (S428) for WB and ELISA. Reacts with human, mouse, and rat STK11. Ideal for LKB1 signaling studies in cancer, metabolism, and cell polarity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-STK11-S428 Rabbit pAb
Remarks/Alias	PJS; LKB1; hLKB1; Phospho-STK11-S428
Species	Human
Gene ID (Human)	6794
Gene ID	6794
Immunogen	A synthetic phosphorylated peptide around S428 of human STK11 (NP_000446.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15831
Protein Weight	49kDa
Shipping	Ice bag