

Phospho-Syk-Y323 Rabbit Pab For Wb, If/Icc - P43405

Item Number: CM0001184



Introduction

Phospho-Syk-Y323 Rabbit Polyclonal Antibody detects human Syk phosphorylated at tyrosine 323. Suitable for Western blot (WB) and immunofluorescence (IF/ICC), it targets spleen tyrosine kinase involved in immune signaling and inflammation.

[Learn More](#)

Parameter	Details
Product Name	Phospho-Syk-Y323 Rabbit pAb
Remarks/Alias	IMD82; p72-Syk; Phospho-Syk-Y323
Species	Human
GeneID (human)	6850
GeneID	6850
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC
Cross Reactivity	Human
SWISS	P43405
Protein Weight	72kDa
Shipping	Ice bag