


Phospho-Tak1 (S412) Rabbit Polyclonal Antibody - O43318

Item Number: CM0001701



Introduction

Phospho-TAK1 (S412) rabbit polyclonal antibody validated for WB and ELISA in human, mouse, rat samples. Targets a key MAP3K involved in TGF- β , IL-1, and TLR signaling pathways, essential for NF- κ B and JNK activation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TAK1-S412 Rabbit pAb
Remarks/Alias	CSCF; FMD2; TAK1; MEKK7; TGF1a; Phospho-TAK1-S412
Species	Human
GeneID (Human)	6885
GeneID	6885
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43318
Protein Weight	67kDa
Shipping	Ice bag