

Phospho-Tau (S214) Rabbit Mab For Wb, Ihc-P, Elisa - P10636

Item Number: CM0000904



Introduction

Phospho-Tau (S214) Rabbit mAb for detection of phosphorylated Tau in human and mouse samples by WB, IHC-P, ELISA. Ideal for neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau (S214) Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau (S214)
Species	Human
GeneID (Human)	4137
Immunogen	Synthetic peptide: A phospho synthetic peptide corresponding to residues surrounding S214 of Human Tau.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag