

Phospho-Tau-S262 Rabbit Pab - P10636

Item Number: CM0000530



Introduction

Phospho-Tau-S262 Rabbit pAb detects Tau phosphorylated at Ser262. Cross-reacts with human, mouse, and rat. Validated for WB, IF/ICC, IF-P, and ELISA. Use in neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-S262 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-S262
Species	Human
GeneID (human)	4137
GeneID	4137
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S262 of Tau (NP_058518.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag