

Phospho-Tau-S404 Rabbit Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P10636

Item Number: CM0000554



Introduction

Rabbit polyclonal antibody against phosphorylated Tau (Ser404) validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with Human, Mouse, Rat. Suitable for neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-S404 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-S404
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	A synthetic phosphorylated peptide around S404 of human Tau (NP_005901.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79 kDa
Shipping	Ice bag