

Phospho-Tau (Ser198/Ser199/Ser202) Rabbit Polyclonal Antibody For Wb And Elisa - P10636

Item Number: CM0001557



Introduction

Rabbit polyclonal antibody targeting human Tau phosphorylated at Ser198/Ser199/Ser202. Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for studying tauopathies and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-Ser198/Ser199/Ser202 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-Ser198/Ser199/Ser202
Species	Human
Gene ID (Human)	4137
Gene ID	4137
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag