

Phospho-Tau-Ser404 Rabbit Pab - P10636

Item Number: CM0001907



Introduction

Phospho-Tau-Ser404 Rabbit polyclonal antibody validated for WB, IF-P, ELISA. Detects human, mouse, rat tau phosphorylated at Ser404. Suitable for Alzheimer's disease and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-Ser404 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-Ser404
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S404 of human Tau (NP_005901.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag