

Phospho-Tau-T181 Rabbit Mab For Wb, Elisa - P10636

Item Number: CM0001891



Introduction

Phospho-Tau-T181 Rabbit monoclonal antibody specifically detects human Tau phosphorylated at T181. Ideal for Western Blot and ELISA applications in neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T181 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T181
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T181 of human Phospho-Tau-T181 (NP_001190181.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag