

Phospho-Tau-T181 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Ihc-P, Elisa - P10636

Item Number: CM0001151



Introduction

Rabbit monoclonal antibody targeting phosphorylated Tau at T181 (SWISS P10636). Validated for WB, IP, IF-P, IHC-P, ELISA in human, mouse, rat. Ideal for Alzheimer's and neuropathology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T181 Rabbit mAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T181
Species	Human
GeneID (human)	4137
GeneID	4137
Immunogen	A synthetic phosphorylated peptide around T181 of human Tau
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag