

Phospho-Tau-T181 Rabbit Pab For Wb And Elisa - P10636-5

Item Number: CM0000484

Introduction

Phospho-Tau-T181 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human and mouse phosphorylated Tau at T181, suitable for Alzheimer's disease and tauopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T181 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T181
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10636-5
Protein Weight	32-45kDa/78-80kDa
Shipping	Ice bag