

Phospho-Tau-T181 Rabbit Polyclonal Antibody (Wb, If, Elisa) - P10636

Item Number: CM0001056



Introduction

High-specificity rabbit polyclonal antibody against Tau phosphorylated at Thr181. Validated for WB, IF/ICC, ELISA; reacts with human, mouse, rat. Essential for Alzheimer's and tauopathy IVD research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T181 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T181
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T498 of human Tau (NP_058519.3?)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag