

Phospho-Tau-T231 Rabbit Pab For Wb, Elisa - P10636

Item Number: CM0001205



Introduction

Phospho-Tau-T231 rabbit polyclonal antibody validated for WB and ELISA. Recognizes mouse and rat species. Targets the Alzheimer's-related phospho-epitope of human Tau protein (79kDa). Ideal for neuroscience and neurodegenerative disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T231 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T231
Species	Human
GeneID (Human)	4137
GeneID	4137
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag