

Phospho-Tau-T231 Rabbit Pab - P10636

Item Number: CM0000701



Introduction

Phospho-Tau-T231 Rabbit pAb (P10636) specifically detects phosphorylated Tau at T231 in human and mouse. Suitable for WB and ELISA. Ideal for Alzheimer's disease and tauopathy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tau-T231 Rabbit pAb
Remarks/Alias	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; tau-40; FTDP-17; PPP1R103; Tau-PHF6; Phospho-Tau-T231
Species	Human
Gene ID	4137
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T548 of human Tau (NP_058519.3?)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P10636
Protein Weight	79kDa
Shipping	Ice bag