


## Phospho-Tifa-T9 Rabbit Pab For Wb And Elisa - Q96Cg3

Item Number: CM0001897



### Introduction

Rabbit polyclonal antibody specifically recognizing human TIFA phosphorylated at Thr-9, validated for Western blot and ELISA, detecting the 21 kDa innate immunity adapter protein.

[Learn More](#)

| Parameter        | Value                                                                    |
|------------------|--------------------------------------------------------------------------|
| Product Name     | Phospho-TIFA-T9 Rabbit pAb                                               |
| Remarks/Alias    | T2BP; T6BP; TIFAA; Phospho-TIFA-T9                                       |
| Species          | Human                                                                    |
| Gene ID (Human)  | 92610                                                                    |
| Immunogen        | A synthetic phosphorylated peptide around T9 of human TIFA (NP_443096.1) |
| Source           | Rabbit                                                                   |
| Category         | Phosphorylated Antibodies                                                |
| Application      | WB, ELISA                                                                |
| Cross Reactivity | Human                                                                    |
| SWISS            | Q96CG3                                                                   |
| Protein Weight   | 21 kDa                                                                   |
| Shipping         | Ice bag                                                                  |