

Phospho-Tp73-Y99 Rabbit Polyclonal Antibody For Wb, Elisa - O15350

Item Number: CM0000692



Introduction

Phospho-TP73-Y99 Rabbit pAb for Western blot and ELISA. Detects human phospho-p73 at Tyr99. Ideal for apoptosis, DNA damage, and lung differentiation research. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-TP73-Y99 Rabbit pAb
Remarks/Alias	P73; CILD47; Phospho-TP73-Y99
Species	Human
GeneID (Human)	7161
GeneID	7161
Immunogen	Synthetic peptide, a phospho specific peptide corresponding to residues surrounding Y99 of human TP73
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15350
Protein Weight	70kDa
Shipping	Ice bag