


Phospho-Tyk2 (Y1054/1055) Rabbit Polyclonal Antibody For Wb And Elisa - P29597

Item Number: CM0001799



Introduction

Phospho-Tyk2 (Y1054/1055) rabbit polyclonal antibody, validated for WB and ELISA, detects human, mouse, and rat Tyk2. Ideal for investigating cytokine and interferon signaling, cell growth, and immune response.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tyk2-Y1054/1055 Rabbit pAb
Remarks/Alias	JTK1; IMD35; Phospho-Tyk2-Y1054/1055
Species	Human
Gene ID (Human)	7297
Gene ID	7297
Immunogen	Synthetic phosphorylated peptide around Y1054 & Y1055 of human Tyk2 (NP_003322.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P29597
Protein Weight	134kDa
Shipping	Ice bag