


Phospho-Tyrosine Hydroxylase-S40 Rabbit Polyclonal Antibody For Wb, Elisa - P07101

Item Number: CM0001383



Introduction

Phospho-Tyrosine Hydroxylase-S40 Rabbit pAb detects phosphorylated tyrosine hydroxylase at Ser40. Suitable for WB and ELISA, cross-reacts with rat. Ideal tool for catecholamine and neurological research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Tyrosine Hydroxylase-S40 Rabbit pAb
Remarks/Alias	TYH; DYT14; DYT5b; Phospho-Tyrosine Hydroxylase-S40
Species	Human
GeneID (Human)	7054
GeneID	7054
Immunogen	Synthetic phosphorylated peptide around S40 of human Tyrosine Hydroxylase (NP_000351.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P07101
Protein Weight	45-59 kDa
Shipping	Ice bag