

## Phospho-Ubiquitin (Ser65) Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P62979 / P62987 / P0Cg47 / P0Cg48

Item Number: CM0004137



### Introduction

Phospho-Ubiquitin (Ser65) Rabbit Monoclonal Antibody validated for WB, IF/ICC, and ELISA. Specific for human ubiquitin phosphorylated at serine 65. Ideal for studying ubiquitin signaling and proteasome pathways.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-Ubiquitin-S65 Rabbit mAb                                                                                    |
| Remarks/Alias    | UBC; S27A; eS31; CEP80; UBA80; HEL112; UBCEP1; UBCEP80; L40; CEP52; RPL40; HUBCEP52; HMG20; Phospho-Ubiquitin-Ser65 |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 6233/7311/7314/7316                                                                                                 |
| Gene ID          | 6233/7311/7314/7316                                                                                                 |
| Immunogen        | A synthetic phosphorylated peptide around S65 of human Ubiquitin (NP_002945.1)                                      |
| Source           | Rabbit                                                                                                              |
| Category         | Monoclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human                                                                                                               |
| SWISS            | P62979/P62987/P0CG47/P0CG48                                                                                         |
| Protein Weight   | 18 kDa/15 kDa/26 kDa/77 kDa                                                                                         |
| Shipping         | Ice bag                                                                                                             |