

Phospho-ULK1-S467 Rabbit Pab For Wb, Elisa - 075385

Item Number: CM0001884



Introduction

Phospho-specific rabbit polyclonal antibody targeting ULK1 (S467), a key autophagy kinase. Validated for Western blot (WB) and ELISA with human and mouse reactivity. Ideal for autophagy and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ULK1-S467 Rabbit pAb
Remarks/Alias	ATG1; ATG1A; UNC51; hATG1; Unc51.1; Phospho-ULK1-S467
Species	Human
GeneID (Human)	8408
GeneID	8408
Immunogen	A synthetic phosphorylated peptide around S467 of human ULK1 (NP_003556.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	075385
Protein Weight	113kDa
Shipping	Ice bag