

Phospho-ULK1-S555 Rabbit Pab For Wb, Elisa - 075385

Item Number: CM0001170



Introduction

Phospho-ULK1-S555 Rabbit pAb specifically detects phosphorylated ULK1 at Ser555. Validated for WB and ELISA in human samples. Ideal for autophagy research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ULK1-S555 Rabbit pAb
Remarks/Alias	ATG1; ATG1A; UNC51; hATG1; Unc51.1
Species	Human
GeneID	8408
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S555 of human ULK1 (NP_003556.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075385
Protein Weight	113 kDa
Shipping	Ice bag