

Phospho-ULK1-S638 Rabbit Pab For Wb And Elisa - O75385

Item Number: CM0001341



Introduction

Rabbit polyclonal antibody targeting phospho-ULK1 (Ser638) for autophagy research. Validated for Western blot and ELISA, cross-reacts with human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ULK1-S638 Rabbit pAb
Remarks/Alias	ATG1; ATG1A; UNC51; hATG1; Unc51.1; Phospho-ULK1-S638
Species	Human
GeneID (human)	8408
GeneID	8408
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S638 of human ULK1 (NP_003556.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O75385
Protein Weight	113kDa
Shipping	Ice bag