

Phospho-B-Catenin (S33/S37/T41) Rabbit Polyclonal Antibody For IHC-P, Elisa - P35222

Item Number: CM0001129



Introduction

Rabbit polyclonal antibody targeting phosphorylated β -Catenin at S33, S37, and T41 for IHC-P and ELISA. Reacts with human, mouse, and rat. Suitable for Wnt signaling and cell adhesion studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho- β -Catenin-S33/S37/T41 Rabbit pAb
Remarks/Alias	EV7; CTNNB; MRD19; NEDSDV; armadillo; Phospho- β -Catenin-S33/S37/T41
Species	Human
GeneID (Human)	1499
GeneID	1499
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag