

Phospho-B-Catenin-S552 Rabbit Monoclonal Antibody For Wb, Elisa - P35222

Item Number: CM0001939



Introduction

Rabbit monoclonal antibody against phospho- β -catenin (Ser552) for Western blot and ELISA applications, with cross-reactivity to human and mouse. Ideal for studying Wnt signaling and cell adhesion.

[Learn More](#)

Parameter	Specification
Product Name	Phospho- β -Catenin-S552 Rabbit mAb
Remarks/Alias	EVR7; CTNNB; MRD19; NEDSDV; armadillo; Phospho- β -Catenin-S552
Species	Human
Gene ID (Human)	1499
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around S552 of human CTNNB1 (NP_001895.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag