

Phospho-B-Catenin-S675 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P35222

Item Number: CM0001612



Introduction

Rabbit polyclonal antibody detecting phospho- β -catenin (S675) for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for Wnt signaling and cell adhesion research.

[Learn More](#)

Parameter	Specification
Product Name	Phospho- β -Catenin-S675 Rabbit pAb
Remarks/Alias	EV7; CTNNB; MRD19; NEDSDV; armadillo; Phospho- β -Catenin-S675
Species	Human
GeneID (Human)	1499
GeneID	1499
Immunogen	A synthetic phosphorylated peptide around S675 of human beta Catenin (NP_001091679.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag