

Phospho-β-Catenin-Y333 Rabbit Pab For Wb - P35222

Item Number: CM0000528



Introduction

Rabbit polyclonal antibody specific for phosphorylated β-Catenin (Tyr333), validated for Western blot in human samples. Ideal for Wnt signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-β-Catenin-Y333 Rabbit pAb
Remarks/Alias	EV7; CTNNB; MRD19; NEDSDV; armadillo; Phospho-β-Catenin-Y333
Species	Human
Gene ID (Human)	1499
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P35222
Protein Weight	85kDa
Shipping	Ice bag