

Phospho-Vangl2-T78/S79/S82 Rabbit Mab For Wb And Elisa - Q9UIK5

Item Number: CM0000655



Introduction

Phospho-specific monoclonal antibody targeting VANGL2 at T78/S79/S82, validated for WB and ELISA. Suitable for human, mouse, and rat planar cell polarity and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VANGL2-T78/S79/S82 Rabbit mAb
Remarks/Alias	LPP1; LTAP; STB1; STBM; STBM1; Phospho-VANGL2-T78/S79/S82
Species	Human
Gene ID (Human)	57216
Gene ID	57216
Immunogen	Synthetic Peptide A phospho specific peptide corresponding to residues surrounding T78/S79/S82 of human VANGL2
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULK5
Protein Weight	60kDa
Shipping	Ice bag