

Phospho-Vasp-S157 Rabbit Mab For Wb, Elisa - P50552

Item Number: CM0000950



Introduction

Recombinant rabbit monoclonal antibody against phospho-VASP (Ser157). Validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for studying actin dynamics, cell migration, and platelet activation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VASP-S157 Rabbit mAb
Remarks/Alias	VASP; Phospho-VASP-S157
Species	Human
GeneID (Human)	7408
GeneID	7408
Immunogen	A synthetic phosphorylated peptide around S157 of human VASP (NP_003361.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P50552
Protein Weight	40kDa
Shipping	Ice bag