

## Phospho-Vasp-S239 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P50552

Item Number: CM0001906



### Introduction

Phospho-VASP-S239 rabbit monoclonal antibody, validated for WB, IF/ICC, and ELISA. Recognizes human VASP phosphorylated at serine 239. Ideal for cytoskeleton dynamics, platelet activation, and cell migration research.

[Learn More](#)

| Parameter        | Value                                                                                           |
|------------------|-------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-VASP-S239 Rabbit mAb                                                                    |
| Remarks/Alias    | VASP; vasodilator-stimulated phosphoprotein; Phospho-VASP-S239                                  |
| Species          | Human                                                                                           |
| GeneID           | 7408                                                                                            |
| Immunogen        | Synthetic peptide (A synthetic phosphorylated peptide around S239 of human VASP (NP_003361.1).) |
| Source           | Rabbit                                                                                          |
| Category         | Phosphorylated Antibodies                                                                       |
| Application      | WB, IF/ICC, ELISA                                                                               |
| Cross Reactivity | Human                                                                                           |
| SWISS            | P50552                                                                                          |
| Protein Weight   | 40kDa                                                                                           |
| Shipping         | Ice bag                                                                                         |