


Phospho-Vav1-Y174 Rabbit Pab For Wb - P15498

Item Number: CM0001407



Introduction

Phospho-specific rabbit polyclonal antibody targeting human VAV1 at tyrosine 174 (Y174). Ideal for WB in human and rat samples, exploring Vav protein function in hematopoietic signaling and cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VAV1-Y174 Rabbit pAb
Remarks/Alias	VAV; Phospho-VAV1-Y174
Species	Human
Gene ID	7409
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding Y174 of human VAV1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Rat
SWISS	P15498
Protein Weight	98kDa
Shipping	Ice bag