

Phospho-Vegf Receptor 2-Y1175 Rabbit Mab For Wb, Elisa - P35968

Item Number: CM0000973



Introduction

Phospho-specific rabbit monoclonal antibody targeting VEGFR2 (KDR, FLK1, CD309) phosphorylated at Y1175. Validated for Western blot and ELISA. Ideal for studying angiogenesis and vascular signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-VEGF Receptor 2-Y1175 Rabbit mAb
Remarks/Alias	FLK1; CD309; VEGFR; VEGFR2
Species	Human
Gene ID (Human)	3791
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y1175 of human VEGF Receptor 2 (NP_002244.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P35968
Protein Weight	152 kDa
Shipping	Ice bag