

Phospho-Vimentin-S39 Rabbit Pab For Wb, Ihc-P, Elisa - P08670

Item Number: CM0001253



Introduction

Phospho-Vimentin-S39 rabbit polyclonal antibody for WB, IHC-P, ELISA applications. Detects human, mouse, rat vimentin phosphorylated at Ser39. Ideal for cell migration and cytoskeleton research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Vimentin-S39 Rabbit pAb
Remarks/Alias	CTRCT30; HEL113; Vimentin; VIM; vimentin; Phospho-Vimentin-S39
Species	Human
GeneID (Human)	7431
GeneID	7431
Immunogen	A synthetic phosphorylated peptide around S39 of human VIM (NP_003371.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08670
Protein Weight	54 kDa
Shipping	Ice bag