


Phospho-Vimentin-S56 Rabbit Pab For Wb, Elisa - P08670

Item Number: CM0000601



Introduction

Phospho-Vimentin-S56 Rabbit pAb detects Vimentin phosphorylated at Ser56 in human and rat samples via WB and ELISA. Ideal for studying intermediate filament dynamics, cell migration, and EMT.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Vimentin-S56 Rabbit pAb
Remarks/Alias	CTRCT30; HEL113; Vimentin; VIM; vimentin; Phospho-Vimentin-S56
Species	Human
GeneID	7431
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S56 of human VIM (NP_003371.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P08670
Protein Weight	54kDa
Shipping	Ice bag