


Phospho-Yap1-S127 Rabbit Mab - P46937

Item Number: CM0001223



Introduction

Phospho-YAP1-S127 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for studying Hippo pathway, organ size control, and tumor suppression. Approx. 54 kDa.

[Learn More](#)

Parameter	Value
Product Name	phospho-YAP1-S127 Rabbit mAb
Remarks/Alias	YAP; YKI; COB1; YAP2; YAP-1; YAP65; phospho-YAP1-S127
Species	Human
GeneID (human)	10413
GeneID	10413
Immunogen	Synthetic peptide: A synthesized peptide derived from human YAP1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	54kDa
Shipping	Ice bag