


Phospho-Yap1 (S397) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P46937

Item Number: CM0001944



Introduction

Rabbit polyclonal antibody specific for YAP1 phosphorylated at Ser397. Validated for WB, IF/ICC, ELISA with human, mouse, rat reactivity. Essential for Hippo pathway and cancer research.

[Learn More](#)

Parameter	Value
Product Name	phospho-YAP1-S397 Rabbit pAb
Remarks/Alias	YAP; YKI; COB1; YAP2; YAP-1; YAP65; phospho-YAP1-S397
Species	Human
GeneID (human)	10413
GeneID	10413
Immunogen	A synthetic phosphorylated peptide around S397 of human YAP1 (NP_001123617.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P46937
Protein Weight	54kDa
Shipping	Ice bag