

Phospho-Ywhaz (Ser58)/Ywhah (Ser59)/Ywhag (Ser59)/Ywhab (Ser60) Rabbit Monoclonal Antibody For Wb, Elisa - P63104 / Q04917 / P61981 / P31946

Item Number: CM0001935



Introduction

Phospho-specific rabbit monoclonal antibody targeting 14-3-3 isoforms zeta (Ser58), eta (Ser59), gamma (Ser59), and beta (Ser60). Validated for WB and ELISA with human and mouse reactivity. Ideal for studying 14-3-3 signaling and phosphorylation events.

[Learn More](#)

Parameter	Value
Product Name	Phospho-YWHAZ -Ser58+YWHAH -Ser59+YWHAG-Ser59+YWHAB -Ser60 Rabbit mAb
Remarks/Alias	HEL4; YWHAD; KCIP-1; HEL-S-3; POPCHAS; HEL-S-93; 14-3-3-zeta; YWHA1; 14-3-3GAMMA; DEE56; EIEE56; PPP1R170; GW128; HEL-S-1; HS1; KCIP-1; YWHAA
Species	Human
Gene ID (Human)	7534/7533/7532/7529
Gene ID	7534/7533/7532/7529
Immunogen	Synthetic phospho-specific peptide corresponding to residues surrounding Ser58&Ser59&Ser59&Ser60 of human YWHAZ/YWHAH/YWHAG/YWHAB
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P63104/Q04917/P61981/P31946
Protein Weight	28kDa
Shipping	Ice bag