


Phospho-Zap70-Y493 Rabbit Mab - P43403

Item Number: CM0001386



Introduction

Phospho-ZAP70 (Tyr493) Rabbit Monoclonal Antibody validated for IF/ICC and ELISA. Reacts with human ZAP70, a key tyrosine kinase in T-cell receptor signaling and adaptive immunity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-ZAP70-Y493 Rabbit mAb
Remarks/Alias	ZAP70; ADMIO2; IMD48; SRK; STCD; STD; TZK; ZAP-70; tyrosine-protein kinase ZAP-70; Phospho-ZAP70-Y493
Species	Human
Gene ID	7535
Immunogen	A synthetic phosphorylated peptide around Y493 of human ZAP70
Source	Rabbit
Category	Phosphorylated Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag